

OSTEOGRAM® REPORT

CompuMed Labs
1230 Rosecrans Ave, Ste100
Manhattan Beach, CA 90260
(310) 258-5000

Mon. Jul 29 2002 10:25

DR. John Radiologist M.D.
CompuMed4BMS
5777 W. Century Blvd., Ste 1285
Los Angeles CA 90045

LOW BONE DENSITY
T-Score¹ = -2.7


Name: Today Test4BMS
Address:
City:
State: ZIP Code:
ID: 25
S.S.#: 123-12-1234
Age at Menopause:
Birth Date: 7/29/64
Test#: 1028
Comment:

Sex: F
Ethnic: W
Height: 5'8"
Weight: 164
Age: 38
X-ray Date: 7/29/02

Date	BMD ² Index	T-Score	Z-Score
07/29/02	82.1	-2.6	-2.6
07/29/02	81.4	-2.7	-2.7

- (1) T-Score = Standard deviations above or below mean BMD for young normals (20-50 years old)
(2) BMD = Bone Mineral Density

SKELETAL STATUS



Physician Notes:

T-Score Range

-1.0 and above
-1.0 to -2.5
-2.5 and below

Skeletal Status

Normal
Low bone mass (Osteopenia)
Osteoporosis